

产品名称: 14-3-3 epsilon Rabbit Monoclonal Antibody
产品货号: AMRe03986



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 epsilon Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAE
别名 (Alternative Names)	14-3-3 epsilon; 14-3-3E; 143E; KCIP-1; MDCR; Protein kinase C inhibitor protein-1; YWHAE
基因 ID (Gene ID)	7531
蛋白 ID (SwissProt ID)	P62258.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 29 kDa

研究背景 (Background)

产品名称: 14-3-3 epsilon Rabbit Monoclonal Antibody
产品货号: AMRe03986




Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathway. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of 14-3-3 epsilon in 293T, SW480 lysates using 14-3-3 epsilon antibody.

注意事项 (Note)

For research use only .