

产品名称: Sulfotransferase 2A1 Rabbit Monoclonal Antibody  
产品货号: AMRe03998



## 产品概述 (Summary)

|                        |                                                 |
|------------------------|-------------------------------------------------|
| 产品名称 (Production Name) | Sulfotransferase 2A1 Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody          |
| 宿主 (Host)              | Rabbit                                          |
| 应用 (Application)       | WB,IHC,ICC/IF,FC                                |
| 种属反应性 (Reactivity)     | Human, Mouse, Rat                               |

## 产品性能 (Performance)

|                     |                                                                                                                |
|---------------------|----------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                                            |
| 克隆 (Clonality)      | Monoclonal                                                                                                     |
| 形式 (Form)           | Liquid                                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                       |
| 储存溶液 (Buffer)       | Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein. |
| 纯化方式 (Purification) | Affinity Purification                                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                                   |
|------------------------|---------------------------------------------------|
| 基因名 (Gene Name)        | SULT2A1                                           |
| 别名 (Alternative Names) | HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST |
| 基因 ID (Gene ID)        | 6822                                              |
| 蛋白 ID (SwissProt ID)   | Q06520.                                           |

## 产品应用 (Application)

|                          |                                                                |
|--------------------------|----------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:34 kDa;Observed MW: 34 kDa                       |

## 研究背景 (Background)

Catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands.


产品名称: Sulfotransferase 2A1 Rabbit Monoclonal Antibody  
产品货号: AMRe03998



## 研究领域 (Research Area)

Signal Transduction

## 图片 (Image Data)



Western blot analysis of Sulfotransferase 2A1 in Human fetal liver lysates using Sulfotransferase 2A1 antibody.

## 注意事项 (Note)

For research use only .