

产品名称: AGTR2 Receptor Rabbit Monoclonal Antibody
产品货号: AMRe04001



产品概述 (Summary)

产品名称 (Production Name)	AGTR2 Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGTR2 AGTR2; Angiotensin ii receptor type 2; Angiotensin II type-2 receptor; AT2 receptor; AT2; At2r; MRX88; Angiotensin ii at2; Angiotensin receptor 2; Type-2 angiotensin II receptor; ATGR2
别名 (Alternative Names)	
基因 ID (Gene ID)	186
蛋白 ID (SwissProt ID)	P50052.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:41 kDa;Observed MW: 41 kDa

产品名称: AGTR2 Receptor Rabbit Monoclonal Antibody
产品货号: AMRe04001




研究背景 (Background)

Receptor for angiotensin II. Cooperates with MTUS1 to inhibit ERK2 activation and cell proliferation.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of AGTR2 Receptor in HepG2, 293T, HeLa lysates using AGTR2 Receptor antibody.

注意事项 (Note)

For research use only .