

产品名称: p53 DINP1 Rabbit Monoclonal Antibody
产品货号: AMRe04002



产品概述 (Summary)

产品名称 (Production Name)	p53 DINP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TP53INP1 TP53INP1; P53DINP1; SIP; Tumor protein p53-inducible nuclear protein 1;
别名 (Alternative Names)	Stress-induced protein; p53-dependent damage-inducible nuclear protein 1; p53DINP1
基因 ID (Gene ID)	94241
蛋白 ID (SwissProt ID)	Q96A56.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:27 kDa;Observed MW: 27 kDa

产品名称: p53 DINP1 Rabbit Monoclonal Antibody
产品货号: AMRe04002




研究背景 (Background)

In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of p53 DINP1 in 3T3 lysates using p53 DINP1 antibody.

注意事项 (Note)

For research use only .