

产品名称: Plexin A1 Rabbit Monoclonal Antibody
产品货号: AMRe04017



产品概述 (Summary)

产品名称 (Production Name)	Plexin A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PLXNA1
别名 (Alternative Names)	NOV; NOV; Plexin-A1; Plexin1; PlexinA1; PLXN1; Plxna1
基因 ID (Gene ID)	5361
蛋白 ID (SwissProt ID)	Q9UIW2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:211 kDa;Observed MW: 211 kDa

研究背景 (Background)

Plays a role in axon guidance, invasive growth and cell migration.


产品名称: Plexin A1 Rabbit Monoclonal Antibody
产品货号: AMRe04017



研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Plexin A1 in K562 lysates using Plexin A1 antibody.

注意事项 (Note)

For research use only .