

产品名称: PITX3/PTX3 Rabbit Monoclonal Antibody
产品货号: AMRe04018



产品概述 (Summary)

产品名称 (Production Name)	PITX3/PTX3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PITX3
别名 (Alternative Names)	Pitx3; PTX3
基因 ID (Gene ID)	5309
蛋白 ID (SwissProt ID)	O75364.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:32 kDa;Observed MW: 32 kDa

研究背景 (Background)

Transcriptional regulator which is important for the differentiation and maintenance of meso-diencephalic dopaminergic

产品名称: PITX3/PTX3 Rabbit Monoclonal Antibody
产品货号: AMRe04018




(mdDA) neurons during development. In addition to its importance during development, it also has roles in the long-term survival and maintenance of the mdDA neurons.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of PITX3/PTX3 in 293T, U87-MG lysates using PITX3/PTX3 antibody.

注意事项 (Note)

For research use only .