

产品名称: Cytokeratin 6 Rabbit Monoclonal Antibody
产品货号: AMRe04020



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT6A/KRT6B/KRT6C CK 6A; CK 6B; CK 6C; CK 6D; CK 6E; CK-6C; CK-6E; K6a keratin; K6b keratin; K6C;
别名 (Alternative Names)	K6c keratin; K6d keratin; K6e keratin; Keratin; KRT6A; KRT6B; KRT6C; KRT6D; KRT6E
基因 ID (Gene ID)	3853
蛋白 ID (SwissProt ID)	P02538.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:60 kDa;Observed MW: 60 kDa

产品名称: Cytokeratin 6 Rabbit Monoclonal Antibody
产品货号: AMRe04020




研究背景 (Background)

There are at least six isoforms of human type II keratin-6 (K6). There are two types of cytoskeletal and microfibrillar keratin, I (acidic) and II (neutral to basic) (40-55 and 56-70 kDa, respectively).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Cytokeratin 6 in HACAT lysates using Cytokeratin 6 antibody.

注意事项 (Note)

For research use only .