

产品名称: SDHA Rabbit Monoclonal Antibody
产品货号: AMRe04021



产品概述 (Summary)

产品名称 (Production Name)	SDHA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SDHA
别名 (Alternative Names)	SDHA; SDH2; SDHF; Succinate dehydrogenase [ubiquinone] flavoprotein subunit; mitochondrial; Flavoprotein subunit of complex II; Fp
基因 ID (Gene ID)	6389
蛋白 ID (SwissProt ID)	P31040.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:73 kDa;Observed MW: 73 kDa

研究背景 (Background)

产品名称: SDHA Rabbit Monoclonal Antibody
产品货号: AMRe04021




Flavoprotein (FP) subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of SDHA in HeLa lysates using SDHA antibody.

注意事项 (Note)

For research use only .