

产品名称: Annexin 4/ANXA4 Rabbit Monoclonal Antibody
产品货号: AMRe04025



产品概述 (Summary)

产品名称 (Production Name)	Annexin 4/ANXA4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA4 AIV; AnnexinA4; AnnexinIV; ANX 4; ANX A4; ANX4; ANXA4; Chromobindin4;
别名 (Alternative Names)	EndonexinI; LipocortinIV; P32.5; P33/41; PAPII; PIG28; PP4X; Protein II; Xanx-4; ZAP36
基因 ID (Gene ID)	307
蛋白 ID (SwissProt ID)	P09525.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:36 kDa; Observed MW: 36 kDa

产品名称: Annexin 4/ANXA4 Rabbit Monoclonal Antibody
产品货号: AMRe04025




研究背景 (Background)

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Annexin 4/ANXA4 in Human placenta, HepG2 lysates using Annexin 4/ANXA4 antibody.

注意事项 (Note)

For research use only .