

产品名称: Flightless I Rabbit Monoclonal Antibody
产品货号: AMRe04027



产品概述 (Summary)

产品名称 (Production Name)	Flightless I Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FLII
别名 (Alternative Names)	FLI; FLIL; Fli1;Protein flightless-1 homolog
基因 ID (Gene ID)	2314.0
蛋白 ID (SwissProt ID)	Q13045.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:145 kDa;Observed MW: 145 kDa

研究背景 (Background)

May play a role as coactivator in transcriptional activation by hormone-activated nuclear receptors (NR) and acts in

产品名称: Flightless I Rabbit Monoclonal Antibody
产品货号: AMRe04027




cooperation with NCOA2 and CARM1. Involved in estrogen hormone signaling. Involved in early embryonic development (By similarity).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Flightless I in HeLa, A549, A431 lysates using Flightless I antibody.

注意事项 (Note)

For research use only .