

产品名称: PRAME Rabbit Monoclonal Antibody
产品货号: AMRe04035



产品概述 (Summary)

产品名称 (Production Name)	PRAME Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRAME
别名 (Alternative Names)	MAPE; OIP4; CT130; OIP-4
基因 ID (Gene ID)	23532
蛋白 ID (SwissProt ID)	P78395.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa; Observed MW: 58 kDa

研究背景 (Background)

Functions as a transcriptional repressor, inhibiting the signaling of retinoic acid through the retinoic acid receptors RARA,

产品名称: PRAME Rabbit Monoclonal Antibody
产品货号: AMRe04035




RARB and RARG. Prevents retinoic acid-induced cell proliferation arrest, differentiation and apoptosis.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse testis using PRAME antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PRAME in K562 lysates using PRAME antibody.

注意事项 (Note)

For research use only .