

产品名称: TAF15 Rabbit Monoclonal Antibody
产品货号: AMRe04036



产品概述 (Summary)

产品名称 (Production Name)	TAF15 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAF15 TAF15; RBP56; TAF2N; TATA-binding protein-associated factor 2N; 68 kDa
别名 (Alternative Names)	TATA-binding protein-associated factor; TAF(II)68; TAFII68; RNA-binding protein 56
基因 ID (Gene ID)	8148
蛋白 ID (SwissProt ID)	Q92804.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 77 kDa

产品名称: TAF15 Rabbit Monoclonal Antibody
产品货号: AMRe04036




研究背景 (Background)

RNA and ssDNA-binding protein that may play specific roles during transcription initiation at distinct promoters. Can enter the preinitiation complex together with the RNA polymerase II (Pol II).


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of TAF15 in K562, rat Brain, C6, 3T3 lysates using TAF15 antibody.



Immunocytochemistry analysis of TAF15(green) in HeLa using TAF15 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TAF15 antibody. High-pressure and temperature

产品名称: TAF15 Rabbit Monoclonal Antibody
产品货号: AMRe04036



Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .