

产品名称: PDGFR beta Rabbit Monoclonal Antibody
产品货号: AMRe04053



产品概述 (Summary)

产品名称 (Production Name)	PDGFR beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDGFRB PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri
别名 (Alternative Names)	
基因 ID (Gene ID)	5159
蛋白 ID (SwissProt ID)	P09619.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 124 kDa; Observed MW: 190 kDa

产品名称: PDGFR beta Rabbit Monoclonal Antibody
产品货号: AMRe04053




研究背景 (Background)

This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. These growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of PDGFR beta in SH-SY5Y, C6 lysates using PDGFR beta antibody.



Western blot analysis of PDGFR beta in 3T3 lysates using PDGFR beta antibody.

注意事项 (Note)

For research use only .