

产品名称: ATE1 Rabbit Monoclonal Antibody
产品货号: AMRe04055



产品概述 (Summary)

产品名称 (Production Name)	ATE1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATE1
别名 (Alternative Names)	Arginyltransferase 1; R-transferase 1
基因 ID (Gene ID)	11101
蛋白 ID (SwissProt ID)	O95260.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa; Observed MW: 59 kDa

研究背景 (Background)

Involved in the post-translational conjugation of arginine to the N-terminal aspartate or glutamate of a protein. This

产品名称: ATE1 Rabbit Monoclonal Antibody
产品货号: AMRe04055




arginylation is required for degradation of the protein via the ubiquitin pathway. Does not arginylate cysteine residues .


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of ATE1 in C6 lysates using ATE1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using ATE1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .