

产品名称: Catalase Rabbit Monoclonal Antibody
产品货号: AMRe04059



产品概述 (Summary)

产品名称 (Production Name)	Catalase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cat
别名 (Alternative Names)	Cas1; Cs-1; Cas-1; 2210418N07
基因 ID (Gene ID)	12359.0
蛋白 ID (SwissProt ID)	P24270.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa; Observed MW: 60 kDa

研究背景 (Background)

Occurs in almost all aerobically respiring organisms and serves to protect cells from the toxic effects of hydrogen peroxide.

产品名称: Catalase Rabbit Monoclonal Antibody
产品货号: AMRe04059




Promotes growth of cells.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Catalase in K562, rat Brain, C6, 3T3 lysates using Catalase antibody.



Immunocytochemistry analysis of Catalase (green) in Raw264.7 using Catalase antibody, and DAPI (blue)

注意事项 (Note)

For research use only .