

产品名称: Phospho-HSF1 (Ser326) Rabbit Monoclonal Antibody
产品货号: AMRe04083



产品概述 (Summary)

产品名称 (Production Name)	Phospho-HSF1 (Ser326) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSF1
别名 (Alternative Names)	HSF1; HSTF1; Heat shock factor protein 1; HSF 1; Heat shock transcription factor 1; HSTF 1
基因 ID (Gene ID)	3297
蛋白 ID (SwissProt ID)	Q00613.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 82 kDa

产品名称: Phospho-HSF1 (Ser326) Rabbit Monoclonal Antibody
产品货号: AMRe04083




研究背景 (Background)

DNA-binding protein that specifically binds heat shock promoter elements (HSE) and activates transcription. In higher eukaryotes, HSF is unable to bind to the HSE unless the cells are heat shocked.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)




Western blot analysis of Phospho-HSF1 (Ser326) in K562, HeLa lysates using Phospho-HSF1 (Ser326) antibody.



Immunocytochemistry analysis of HSF1 (Phospho-Ser326)(green) in HeLa using HSF1 (Phospho-Ser326) antibody, and DAPI(blue)

产品名称: Phospho-HSF1 (Ser326) Rabbit Monoclonal Antibody
产品货号: AMRe04083



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using HSF1 (Phospho-Ser326) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .