

产品名称: ICAD Rabbit Monoclonal Antibody  
产品货号: AMRe04085

## 产品概述 (Summary)

产品名称 (Production Name)	ICAD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

## 产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

## 免疫原信息 (Immunogen)

基因名 (Gene Name)	DFFA
别名 (Alternative Names)	DFFA; DFF1; DFF45; H13; DNA fragmentation factor subunit alpha; DNA fragmentation factor 45 kDa subunit; DFF-45; Inhibitor of CAD; ICAD
基因 ID (Gene ID)	1676
蛋白 ID (SwissProt ID)	O00273.

## 产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 45 kDa

## 研究背景 (Background)

产品名称: ICAD Rabbit Monoclonal Antibody  
产品货号: AMRe04085




Inhibitor of the caspase-activated DNase (DFF40).

## 研究领域 (Research Area)

Cell Biology

## 图片 (Image Data)



Western blot analysis of ICAD in rat Brain, MCF-7, C6 lysates using ICAD antibody.

## 注意事项 (Note)

For research use only .