

产品名称: MAD2L2 Rabbit Monoclonal Antibody
产品货号: AMRe04088



产品概述 (Summary)

产品名称 (Production Name)	MAD2L2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAD2L2
别名 (Alternative Names)	REV7; FANCV; MAD2B; POLZ2
基因 ID (Gene ID)	10459
蛋白 ID (SwissProt ID)	Q9UI95.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 24 kDa

研究背景 (Background)

Adapter protein able to interact with different proteins and involved in different biological processes. Mediates the

产品名称: MAD2L2 Rabbit Monoclonal Antibody
产品货号: AMRe04088




interaction between the error-prone DNA polymerase zeta catalytic subunit REV3L and the inserter polymerase REV1, thereby mediating the second polymerase switching in translesion DNA synthesis.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Mad2L2/REV7 in Jurkat lysates using MAD2L2 antibody.



Immunocytochemistry analysis of Mad2L2/REV7 (green) in Hela using Mad2L2/REV7 antibody, and DAPI (blue)

注意事项 (Note)

For research use only .