

产品名称: Nicastrin Rabbit Monoclonal Antibody
产品货号: AMRe04092



产品概述 (Summary)

产品名称 (Production Name)	Nicastrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCSTN
别名 (Alternative Names)	APH2; Ncstn; NCT; Nicastrin
基因 ID (Gene ID)	23385
蛋白 ID (SwissProt ID)	Q92542.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 78 kDa; Observed MW: 110-120 kDa

研究背景 (Background)

Essential subunit of the gamma-secretase complex, an endoprotease complex that catalyzes the intramembrane cleavage

产品名称: Nicastrin Rabbit Monoclonal Antibody
产品货号: AMRe04092




of integral membrane proteins such as Notch receptors and APP (beta-amyloid precursor protein). It probably represents a stabilizing cofactor required for the assembly of the gamma-secretase complex.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Nicastrin in K562, rat Brain, 3T3, HeLa lysates using Nicastrin antibody.

注意事项 (Note)

For research use only .