

产品名称: PI3 Kinase p110 beta Rabbit Monoclonal Antibody
产品货号: AMRe04097



产品概述 (Summary)

产品名称 (Production Name)	PI3 Kinase p110 beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PIK3CB
别名 (Alternative Names)	PIK3CB; DKFZp779K1237; MGC133043; PI3K; PI3KCB; PI3Kbeta; PIK3C1; p110-BETA
基因 ID (Gene ID)	5291
蛋白 ID (SwissProt ID)	P42338.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 123 kDa; Observed MW: 110 kDa

研究背景 (Background)

产品名称: PI3 Kinase p110 beta Rabbit Monoclonal Antibody
产品货号: AMRe04097




Phosphoinositide 3-kinase (PI3K) catalyzes the production of phosphatidylinositol-3,4,5-triphosphate by phosphorylating phosphatidylinositol (PI), phosphatidylinositol-4-phosphate (PIP) and phosphatidylinositol-4,5-bisphosphate (PIP₂). Growth factors and hormones trigger this phosphorylation event, which in turn coordinates cell growth, cell cycle entry, cell migration, and cell survival.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of PI3 Kinase p110 beta in HCT116, SH-SY5Y lysates using PI3 Kinase p110 beta antibody.

注意事项 (Note)

For research use only .