

产品名称: Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
产品货号: AMRe04114



产品概述 (Summary)

产品名称 (Production Name)	Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDRG1
别名 (Alternative Names)	Protein NDRG1; Differentiation-related gene 1 protein; DRG-1; RTP; Rit42; NDRG1; CAP43; DRG1; RTP; targ1; TDD5; tdds
基因 ID (Gene ID)	10397
蛋白 ID (SwissProt ID)	Q92597.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43 kDa

产品名称: Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
产品货号: AMRe04114




研究背景 (Background)

May have a growth inhibitory role.


研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Immunocytochemistry analysis of Phospho-NDRG1 (Ser330) (green) in 293T using Phospho-NDRG1 (Ser330) antibody, and DAPI(blue).




Western blot analysis of Phospho-NDRG1 (Ser330) in mouse kidney lysates using Phospho-NDRG1 (Ser330) antibody.



产品名称: Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
产品货号: AMRe04114



Western blot analysis of Phospho-NDRG1(Ser330) in Hela, A549, MCF-7, 239T, C2C12 lysates using Phospho-NDRG1(Ser330) antibody



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using NDRG1 (Phospho-Ser330) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .