

**产品名称: Annexin 4 Rabbit Monoclonal Antibody**  
**产品货号: AMRe04116**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Annexin 4 Rabbit Monoclonal Antibody   |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF                          |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat,Hamster                |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                                                                                        |
|------------------------|----------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | ANXA4<br>AIV; AnnexinA4; AnnexinIV; ANX 4; ANX A4; ANX4; ANXA4; Chromobindin4;         |
| 别名 (Alternative Names) | EndonexinI; LipocortinIV; P32.5; P33/41; PAPII; PIG28; PP4X; Protein II; Xanx-4; ZAP36 |
| 基因 ID (Gene ID)        | 307                                                                                    |
| 蛋白 ID (SwissProt ID)   | P09525.                                                                                |

## 产品应用 (Application)

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 36 kDa; Observed MW: 36 kDa       |

产品名称: Annexin 4 Rabbit Monoclonal Antibody  
产品货号: AMRe04116




## 研究背景 (Background)

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.


## 研究领域 (Research Area)

Cardiovascular


## 图片 (Image Data)



Immunocytochemistry analysis of Annexin 4 (green) in Hela using Annexin 4 antibody, and DAPI (blue).




Western blot analysis of Annexin4/ANXA4 in mouse kidney, mouse lung lysates using Annexin4/ANXA4 antibody.



Western blot analysis of Annexin A4 in Hela, CHO-K1, C6 lysates using Annexin A4 antibody

**产品名称: Annexin 4 Rabbit Monoclonal Antibody**  
**产品货号: AMRe04116**



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Annexin4/ANXA4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

#### **注意事项 (Note)**

For research use only .