

产品名称: Olig1 Rabbit Monoclonal Antibody
产品货号: AMRe04130



产品概述 (Summary)

产品名称 (Production Name)	Olig1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OLIG1
别名 (Alternative Names)	BHLH B6; BHLHB 6; bHLHb6; bHLHe21; Olig1; Oligo1
基因 ID (Gene ID)	116448
蛋白 ID (SwissProt ID)	Q8TAK6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 35 kDa

研究背景 (Background)

Promotes formation and maturation of oligodendrocytes, especially within the brain. Cooperates with OLIG2 to establish

产品名称: Olig1 Rabbit Monoclonal Antibody
产品货号: AMRe04130




the pMN domain of the embryonic neural tube.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Olig1 in 3T3 lysates using Olig1 antibody.



Western blot analysis of Olig1 in A549, HL-60, U2OS, C6 lysates using Olig1 antibody.

注意事项 (Note)

For research use only .