

产品名称: PAR4 Rabbit Monoclonal Antibody
产品货号: AMRe04131



产品概述 (Summary)

产品名称 (Production Name)	PAR4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAWR
别名 (Alternative Names)	PAWR; PAR4; PRKC apoptosis WT1 regulator protein; Prostate apoptosis response 4 protein; Par-4
基因 ID (Gene ID)	5074
蛋白 ID (SwissProt ID)	Q96IZ0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 41 kDa

研究背景 (Background)

产品名称: PAR4 Rabbit Monoclonal Antibody
产品货号: AMRe04131




Pro-apoptotic protein capable of selectively inducing apoptosis in cancer cells, sensitizing the cells to diverse apoptotic stimuli and causing regression of tumors in animal models. Induces apoptosis in certain cancer cells by activation of the Fas prodeath pathway and coparallel inhibition of NF-kappa-B transcriptional activity.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of PAR4 in 3T3 lysates using PAR4 antibody.



Western blot analysis of PAR4 in HeLa, A549, HL-60, U251, U87-MG lysates using PAR4 antibody.

注意事项 (Note)

For research use only .