

产品名称: MLKL Rabbit Monoclonal Antibody
产品货号: AMRe04148



产品概述 (Summary)

产品名称 (Production Name)	MLKL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLKL
别名 (Alternative Names)	MLKL
基因 ID (Gene ID)	197259
蛋白 ID (SwissProt ID)	Q8NB16.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54 kDa

研究背景 (Background)

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process. Activated following

产品名称: MLKL Rabbit Monoclonal Antibody
产品货号: AMRe04148




phosphorylation by RIPK3, leading to homotrimerization, localization to the plasma membrane and execution of programmed necrosis characterized by calcium influx and plasma membrane damage. Does not have protein kinase activity.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of MLKL in HUVEC, RAW264.7 lysates using MLKL antibody.

注意事项 (Note)

For research use only .