

产品名称: Acetyl Coenzyme A Carboxylase (6N12) Rabbit Monoclonal Antibody
产品货号: AMRe04158



产品概述 (Summary)

产品名称 (Production Name)	Acetyl Coenzyme A Carboxylase (6N12) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% protective protein
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACACB
别名 (Alternative Names)	
基因 ID (Gene ID)	32.0
蛋白 ID (SwissProt ID)	O00763.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:100-1:500
蛋白分子量 (Molecular Weight)	277 kDa


研究背景 (Background)

产品名称: Acetyl Coenzyme A Carboxylase (6N12) Rabbit
Monoclonal Antibody
产品货号: AMRe04158



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Acetyl Coenzyme A Carboxylase in HeLa, HCT116 lysates using Acetyl Coenzyme A Carboxylase antibody.

注意事项 (Note)

For research use only .