

产品名称: Acetyl-Histone H3 (Lys27) (17F16) Rabbit Monoclonal Antibody
产品货号: AMRe04172



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-Histone H3 (Lys27) (17F16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H3C1
别名 (Alternative Names)	
基因 ID (Gene ID)	8350.0
蛋白 ID (SwissProt ID)	P68431.A synthetic acetyl-peptide corresponding to residues surrounding Lys27 of human Histone H3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:50000,IHC 1:1000-1:10000,ICC/IF 1:500-1:5000,FC 1:200-1:1000
蛋白分子量 (Molecular Weight)	15kDa

产品名称: Acetyl-Histone H3 (Lys27) (17F16) Rabbit
Monoclonal Antibody
产品货号: AMRe04172




研究背景 (Background)

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of Histone H3 (acetyl-Lys27) in HeLa,A549,HL-60 cell lysates using Histone H3 (acetyl-Lys27) antibody(1:1000 diluted).

注意事项 (Note)

For research use only .