

产品名称: Phospho-B Raf (T401) (6C16) Rabbit Monoclonal Antibody
产品货号: AMRe05857



产品概述 (Summary)

产品名称 (Production Name)	Phospho-B Raf (T401) (6C16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BRAF
别名 (Alternative Names)	B-RAF; B-RAF proto-oncogene serine/threonine-protein kinase; BRAF1; RAFB1; RMIL; RMIL serine/threonine-protein kinase; c-RMIL; kinase B-Raf;
基因 ID (Gene ID)	673.0
蛋白 ID (SwissProt ID)	P15056.A synthetic phosphopeptide corresponding to residues surrounding Thr401 of human B Raf

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:20-1:50, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	84kDa

产品名称: Phospho-B Raf (T401) (6C16) Rabbit Monoclonal Antibody
产品货号: AMRe05857


研究背景 (Background)

Protein kinase involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron. Phosphorylates MAP2K1, and thereby contributes to the MAP kinase signal transduction pathway. Protein kinase involved in the transduction of mitogenic signals from the cell membrane to the nucleus (Probable). Phosphorylates MAP2K1, and thereby activates the MAP kinase signal transduction pathway (PubMed:21441910, PubMed:29433126). May play a role in the postsynaptic responses of hippocampal neurons (PubMed:1508179).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of Phospho-BRAF (Thr401) in K562,C6,3T3,HeLa cell lysates using Phospho-BRAF (Thr401) antibody(1:1000 diluted).

注意事项 (Note)

For research use only .