

**产品名称: Phospho-CDK1/2 (Thr14) (5Z8) Rabbit
Monoclonal Antibody**
产品货号: AMRe05871



产品概述 (Summary)

产品名称 (Production Name)	Phospho-CDK1/2 (Thr14) (5Z8) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDK1/CDK2/
别名 (Alternative Names)	
基因 ID (Gene ID)	983/1017
蛋白 ID (SwissProt ID)	P06493/P24941.A synthetic phosphopeptide corresponding to residues surrounding Thr14 of human Cdk1/2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,IP 1:20-1:50,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	34kDa

产品名称: Phospho-CDK1/2 (Thr14) (5Z8) Rabbit
Monoclonal Antibody
产品货号: AMRe05871




研究背景 (Background)

研究领域 (Research Area)

Cell Biology; Cell Cycle; Kinases/Phosphatases; Cdks; Neuroscience; Neurology process; Neural Signal Transduction; Neurodegenerative disease; Epigenetics and Nuclear Signaling; Cancer

图片 (Image Data)



Western blot detection of Cdk1/2 (phospho-Thr14) in HepG2,U2OS,U251 using Cdk1/2 (phospho-Thr14) antibody(1:1000 diluted)

注意事项 (Note)

For research use only .