

产品名称: Phospho-ErbB2(Y1221 + Y1222) (11Z13) Rabbit
Monoclonal Antibody
产品货号: AMRe05899

产品概述 (Summary)

产品名称 (Production Name)	Phospho-ErbB2(Y1221 + Y1222) (11Z13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ERBB2 CD340; CerbB2; Erb b2 receptor tyrosine kinase 2; ERBB2; HER2; Herstatin; Human epidermal growth factor receptor 2; MLN19; NEU; NGL; Proto-oncogene Neu; Receptor tyrosine-protein kinase erbB-2; Tyrosine kinase type cell surface receptor HER2; V erb b2 avian erythroblastic leukemia viral oncogene homolog 2; V erb b2 avian erythroblastic leukemia viral oncoprotein 2;
别名 (Alternative Names)	
基因 ID (Gene ID)	2064.0
蛋白 ID (SwissProt ID)	P04626.

产品名称: Phospho-ErbB2(Y1221 + Y1222) (11Z13) Rabbit
Monoclonal Antibody
产品货号: AMRe05899

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	138kDa


研究背景 (Background)

Protein tyrosine kinase that is part of several cell surface receptor complexes, but that apparently needs a coreceptor for ligand binding. Essential component of a neuregulin-receptor complex, although neuregulins do not interact with it alone. GP30 is a potential ligand for this receptor. Protein tyrosine kinase that is part of several cell surface receptor complexes, but that apparently needs a coreceptor for ligand binding. Essential component of a neuregulin-receptor complex, although neuregulins do not interact with it alone. GP30 is a potential ligand for this receptor. Regulates outgrowth and stabilization of peripheral microtubules (MTs). Upon ERBB2 activation, the MEMO1-RHOA-DIAPH1 signaling pathway elicits the phosphorylation and thus the inhibition of GSK3B at cell membrane. This prevents the phosphorylation of APC and CLASP2, allowing its association with the cell membrane. In turn, membrane-bound APC allows the localization of MACF1 to the cell membrane, which is required for microtubule capture and stabilization.

研究领域 (Research Area)

Cancer

图片 (Image Data)



Western blot analysis of Phospho-ErbB2(Y1221 + Y1222) expression in SKBR3 cell lysate.

注意事项 (Note)

For research use only .