

**产品名称: Phospho-Glucocorticoid Receptor (Ser226)
(12W11) Rabbit Monoclonal Antibody**
产品货号: AMRe05911



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Glucocorticoid Receptor (Ser226) (12W11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% protective protein
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NR3C1
别名 (Alternative Names)	
基因 ID (Gene ID)	2908.0
蛋白 ID (SwissProt ID)	P04150.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,IF-P 1:200-1:1000
蛋白分子量 (Molecular Weight)	86 kDa

产品名称: Phospho-Glucocorticoid Receptor (Ser226)
(12W11) Rabbit Monoclonal Antibody
产品货号: AMRe05911




研究背景 (Background)


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Phospho-Glucocorticoid Receptor (Ser226) in HeLa, A549, C6 lysates using Phospho-Glucocorticoid Receptor (Ser226) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Glucocorticoid Receptor (Phospho-Ser226) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .