

产品名称: Phospho-Histone H3 (Ser10/Thr11) (17B15)
Rabbit Monoclonal Antibody
产品货号: AMRe05918



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Histone H3 (Ser10/Thr11) (17B15) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H3C1
别名 (Alternative Names)	
基因 ID (Gene ID)	8350.0
蛋白 ID (SwissProt ID)	P68431.A synthetic phosphopeptide corresponding to residues surrounding Ser10/Thr11 of human Histone H3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:200-1:1000,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	15kDa

产品名称: Phospho-Histone H3 (Ser10/Thr11) (17B15)
Rabbit Monoclonal Antibody
产品货号: AMRe05918




研究背景 (Background)

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of Histone H3 (phospho-Ser10/Thr11) in HeLa,A549,HL-60,C6 using Histone H3 (phospho-Ser10/Thr11) antibody(1:1000 diluted)

注意事项 (Note)

For research use only .