

产品名称: Phospho-IRF3 (Ser386) (9P16) Rabbit Monoclonal Antibody
产品货号: AMRe05928



产品概述 (Summary)

产品名称 (Production Name)	Phospho-IRF3 (Ser386) (9P16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRF3
别名 (Alternative Names)	
基因 ID (Gene ID)	3661.0
蛋白 ID (SwissProt ID)	Q14653.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	47 kDa

研究背景 (Background)


产品名称: Phospho-IRF3 (Ser386) (9P16) Rabbit Monoclonal Antibody
产品货号: AMRe05928



研究领域 (Research Area)

Immunology; Innate Immunity; Cytokines; Interferons; Epigenetics and Nuclear Signaling; Transcription; Other factors; Microbiology; Organism; Virus; RNA Virus; ssRNA positive strand virus; SARS Coronavirus; TLR Signaling

图片 (Image Data)



Western blot analysis of Phospho-IRF3 (S386) in MCF-7 lysates using Phospho-IRF3 (Ser386) antibody.

注意事项 (Note)

For research use only .