

产品名称: Phospho-Jun/JunD (Ser73/Ser100) (12F13) Rabbit
Monoclonal Antibody
产品货号: AMRe05932



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Jun/JunD (Ser73/Ser100) (12F13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	JUN/JUND
别名 (Alternative Names)	
基因 ID (Gene ID)	3725/3727
蛋白 ID (SwissProt ID)	P05412/P17535.A synthetic phosphopeptide corresponding to residues surrounding Ser73 of human Jun

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,IP 1:50-1:200,IF-P 1:200-1:1000
蛋白分子量 (Molecular Weight)	35,36kDa

产品名称: Phospho-Jun/JunD (Ser73/Ser100) (12F13) Rabbit
Monoclonal Antibody
产品货号: AMRe05932




研究背景 (Background)

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of Jun (phospho-Ser73)/JunD (phospho-Ser100) in A549,HL-60,U251,U87-MG cell lysates using Jun (phospho-Ser73)/JunD (phospho-Ser100) antibody(1:1000 diluted).

注意事项 (Note)

For research use only .