

**产品名称: Phospho-Paxillin (Y118) (19J11) Rabbit
Monoclonal Antibody**
产品货号: AMRe05963



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Paxillin (Y118) (19J11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PXN
别名 (Alternative Names)	Paired box protein Pax 1; PAX 1; PAX1; Paxillin alpha; Paxillin; PXN; PXN protein;
基因 ID (Gene ID)	5829.0
蛋白 ID (SwissProt ID)	P49023.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	65kDa

产品名称: Phospho-Paxillin (Y118) (19J11) Rabbit
Monoclonal Antibody
产品货号: AMRe05963




研究背景 (Background)

Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Phospho-Paxillin (Y118) expression in (1) HeLa cell lysate; (2) HeLa cell treated with pervanadate.

注意事项 (Note)

For research use only .