

产品名称: 14-3-3 epsilon (1W4) Rabbit Monoclonal Antibody
产品货号: AMRe06274



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 epsilon (1W4) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAE
别名 (Alternative Names)	14-3-3E; 143E; KCIP-1; MDCR; Protein kinase C inhibitor protein-1; YWHAE;
基因 ID (Gene ID)	7531.0
蛋白 ID (SwissProt ID)	P62258.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, FC 1:100-1:1000
蛋白分子量 (Molecular Weight)	29kDa

产品名称: 14-3-3 epsilon (1W4) Rabbit Monoclonal Antibody
产品货号: AMRe06274




研究背景 (Background)

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathway. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner. Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner (By similarity). Positively regulates phosphorylated protein HSF1 nuclear export to the cytoplasm (PubMed:12917326).

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of 14-3-3 epsilon expression in 293T cell lysate.

注意事项 (Note)

For research use only .