

产品名称: 14-3-3 gamma (13M6) Rabbit Monoclonal Antibody
产品货号: AMRe06275



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 gamma (13M6) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAG
别名 (Alternative Names)	1433G, 143G, KCIP-1, Protein kinase C inhibitor protein-1, YWHAG; gamma polypeptide;
基因 ID (Gene ID)	7532.0
蛋白 ID (SwissProt ID)	P61981.A synthetic peptide of human 14-3-3 gamma/YWHAG

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	28kDa

产品名称: 14-3-3 gamma (13M6) Rabbit Monoclonal Antibody
产品货号: AMRe06275




研究背景 (Background)

Induce target protein conformational changes that modify target protein function. Distinct temporal and spatial expression patterns of 14-3-3 isoforms have been observed in development and in acute response to extracellular signals and drugs, suggesting that 14-3-3 isoforms may perform different functions despite their sequence similarities. Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of 14-3-3 gamma/YWHAG in HL-60,U251 cell lysates using 14-3-3 gamma/YWHAG antibody(1:1000 diluted).

注意事项 (Note)

For research use only .