

产品名称: 41701(11H8) Rabbit Monoclonal Antibody
产品货号: AMRe06324



产品概述 (Summary)

产品名称 (Production Name)	41701(11H8) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	14 3 3 protein beta/alpha; 14 3 3 protein zeta; 14 3 3 zeta; 14-3-3
别名 (Alternative Names)	protein/cytosolic phospholipase A2; GW128; HS1; KCIP 1; YWHAB; YWHAD; YWHAZ;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P31946/P63104.A synthetic peptide of human '14-3-3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	28kDa

产品名称: 41701(11H8) Rabbit Monoclonal Antibody
产品货号: AMRe06324




研究背景 (Background)

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot detection of 41701 in Jurkat,Rat Brain,C6,CHO-K1 cell lysates using 14-3-3 antibody(1:500 diluted).

注意事项 (Note)

For research use only .