

产品名称: 45567(3A2) Rabbit Monoclonal Antibody
产品货号: AMRe06328



产品概述 (Summary)

产品名称 (Production Name)	45567(3A2) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POU2F2
别名 (Alternative Names)	Oct 2; Octamer binding transcription factor 2; OTF2; POU domain; POU2F2;
基因 ID (Gene ID)	5452.0
蛋白 ID (SwissProt ID)	P09086.A synthetic peptide of human Oct-2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:200-1:500,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	51kDa

研究背景 (Background)

Transcription factor that specifically binds to the octamer motif (5'-ATTGTCAT-3'). Regulates transcription in a number of


产品名称: 45567(3A2) Rabbit Monoclonal Antibody
产品货号: AMRe06328



tissues in addition to activating immunoglobulin gene expression. Modulates transcription transactivation by NR3C1, AR and PGR. Transcription factor that specifically binds to the octamer motif (5'-ATTGTCAT-3') (PubMed:2904654, PubMed:7859290). Regulates IL6 expression in B cells with POU2AF1 (By similarity). Regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression (PubMed:2901913, PubMed:2904654). Modulates transcription transactivation by NR3C1, AR and PGR (PubMed:10480874).

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of 44106 in Rat Brain lysates using Oct2 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .