

产品名称: 45569(2Y12) Rabbit Monoclonal Antibody
产品货号: AMRe06329



产品概述 (Summary)

产品名称 (Production Name)	45569(2Y12) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POU5F1
别名 (Alternative Names)	Octamer binding protein 3; Octamer binding protein 4; OCT3; OCT4; OTF3; OTF4; OTF-3; Oct-3; Oct-4; POU5F1;
基因 ID (Gene ID)	5460.0
蛋白 ID (SwissProt ID)	Q01860.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:50-1:200,ICC/IF 1:100-1:200,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	39kDa

产品名称: 45569(2Y12) Rabbit Monoclonal Antibody
产品货号: AMRe06329




研究背景 (Background)

Oct-4 is a transcription factor highly expressed in undifferentiated embryonic stem cells and embryonic germ cells. Forms a trimeric complex with SOX2 on DNA and controls the expression of a number of genes involved in embryonic development such as YES1, FGF4, UTF1 and ZFP206. Critical for early embryogenesis and for embryonic stem cell pluripotency.

Transcription factor that binds to the octamer motif (5'- ATTTGCAT-3'). Forms a trimeric complex with SOX2 or SOX15 on DNA and controls the expression of a number of genes involved in embryonic development such as YES1, FGF4, UTF1 and ZFP206. Critical for early embryogenesis and for embryonic stem cell pluripotency.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Oct4 expression in NCCIT cell lysate.

注意事项 (Note)

For research use only .