

产品名称: Aconitase 2 (15R17) Rabbit Monoclonal Antibody
产品货号: AMRe06512



产品概述 (Summary)

产品名称 (Production Name)	Aconitase 2 (15R17) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACO2
别名 (Alternative Names)	ACO2; Aconitase2; Aconitate hydratase; ACONM; Citrate hydro lyase; ICRD;
基因 ID (Gene ID)	50.0
蛋白 ID (SwissProt ID)	Q99798.A synthetic peptide of human Aconitase 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	85kDa

研究背景 (Background)

Catalyzes the isomerization of citrate to isocitrate via cis-aconitate. Catalyzes the isomerization of citrate to isocitrate via

产品名称: Aconitase 2 (15R17) Rabbit Monoclonal Antibody
产品货号: AMRe06512




cis- aconitate.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot detection of Aconitase 2 in HL-60,C6 using Aconitase 2 antibody(1:1000 diluted)

注意事项 (Note)

For research use only .