

产品名称: ACTR1B (17O10) Rabbit Monoclonal Antibody
产品货号: AMRe06561



产品概述 (Summary)

产品名称 (Production Name)	ACTR1B (17O10) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTR1B
别名 (Alternative Names)	ACTR1B; ARP1B; Beta centractin; CTRN2; PC3;
基因 ID (Gene ID)	10120.0
蛋白 ID (SwissProt ID)	P42025.A synthetic peptide of human ACTR1B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:100-1:200,FC 1:200-1:500,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	42kDa

研究背景 (Background)


产品名称: ACTR1B (17O10) Rabbit Monoclonal Antibody
产品货号: AMRe06561



Component of a multi-subunit complex involved in microtubule based vesicle motility. It is associated with the centrosome.
Component of a multi-subunit complex involved in microtubule based vesicle motility. It is associated with the centrosome.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of ACTR1B in Hela,A549,HL-60,U251,U87-MG cell lysates using ACTR1B antibody(1:1000 diluted).

注意事项 (Note)

For research use only .