

产品名称: AFP (1J18) Rabbit Monoclonal Antibody
产品货号: AMRe06665



产品概述 (Summary)

产品名称 (Production Name)	AFP (1J18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AFP
别名 (Alternative Names)	Alpha fetoglobulin; Alpha fetoprotein precursor; FETA; HPAFP;
基因 ID (Gene ID)	174.0
蛋白 ID (SwissProt ID)	P02771.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,IP 1:20-1:50,IF-P 1:50-1:200
蛋白分子量 (Molecular Weight)	69kDa

研究背景 (Background)

产品名称: AFP (1J18) Rabbit Monoclonal Antibody
产品货号: AMRe06665




Alpha-fetoprotein (AFP) is a 65 kDa glycoprotein found in the mammalian fetal liver, yolk sac, and GI tract. In normal adult cells, p53 binds to the repressor region of the AFP gene, thereby blocking transcription. Mutations in both p53 and β -catenin are associated with aberrant expression of AFP. Research studies have shown that elevated serum AFP levels are predictive of hepatocellular carcinoma .

Binds copper, nickel, and fatty acids as well as, and bilirubin less well than, serum albumin. Only a small percentage (less than 2%) of the human AFP shows estrogen-binding properties.

研究领域 (Research Area)

Cancer

图片 (Image Data)



Western blot analysis of AFP expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .