

产品名称: AHA1 (3L10) Rabbit Monoclonal Antibody
产品货号: AMRe06690



产品概述 (Summary)

产品名称 (Production Name)	AHA1 (3L10) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AHSA1
别名 (Alternative Names)	AHA 1; AHA1; AHSA 1; Ahsa1; C14orf3; HSPC322; p38;
基因 ID (Gene ID)	10598.0
蛋白 ID (SwissProt ID)	O95433.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	38kDa

研究背景 (Background)

产品名称: AHA1 (3L10) Rabbit Monoclonal Antibody
产品货号: AMRe06690




Cochaperone that stimulates HSP90 ATPase activity (By similarity). May affect a step in the endoplasmic reticulum to Golgi trafficking. Acts as a co-chaperone of HSP90AA1 (PubMed:29127155). Activates the ATPase activity of HSP90AA1 leading to increase in its chaperone activity (PubMed:29127155). Competes with the inhibitory co- chaperone FNIP1 for binding to HSP90AA1, thereby providing a reciprocal regulatory mechanism for chaperoning of client proteins (PubMed:27353360). Competes with the inhibitory co-chaperone TSC1 for binding to HSP90AA1, thereby providing a reciprocal regulatory mechanism for chaperoning of client proteins (PubMed:29127155).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of AHA1 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .