

产品名称: alpha Tubulin (4D19) Rabbit Monoclonal Antibody
产品货号: AMRe06803



产品概述 (Summary)

产品名称 (Production Name)	alpha Tubulin (4D19) Rabbit Monoclonal Antibody
描述 (Description)	Rabbit Monoclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUBA4A
别名 (Alternative Names)	
基因 ID (Gene ID)	7277.0
蛋白 ID (SwissProt ID)	P68366.A synthetic peptide of human alpha Tubulin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB: 1:1000-1:10000
蛋白分子量 (Molecular Weight)	50kDa


研究背景 (Background)

产品名称: alpha Tubulin (4D19) Rabbit Monoclonal Antibody
产品货号: AMRe06803



研究领域 (Research Area)

图片 (Image Data)



Western blot detection of alpha Tubulin in HeLa,CHO-K1,Rat Brain lysates using alpha Tubulin antibody(1:500 diluted).

注意事项 (Note)

For research use only .