

产品名称: Angiotensinogen (2C9) Rabbit Monoclonal Antibody
产品货号: AMRe06891



产品概述 (Summary)

产品名称 (Production Name)	Angiotensinogen (2C9) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGT
别名 (Alternative Names)	AGT; Ang; Ang I; Ang II; Ang III; Angiotensinogen; ANHU; ANRT; AT-2; AT-II; PAT; Pre angiotensinogen; Serpin A8; SERPINA8;
基因 ID (Gene ID)	183.0
蛋白 ID (SwissProt ID)	P01019.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000
蛋白分子量 (Molecular Weight)	53kDa

产品名称: Angiotensinogen (2C9) Rabbit Monoclonal Antibody
产品货号: AMRe06891




研究背景 (Background)

Essential component of the renin-angiotensin system (RAS), a potent regulator of blood pressure, body fluid and electrolyte homeostasis. In response to lowered blood pressure, the enzyme renin cleaves angiotensinogen to produce angiotensin-1 (angiotensin 1-10). Essential component of the renin-angiotensin system (RAS), a potent regulator of blood pressure, body fluid and electrolyte homeostasis. [Angiotensin-3]: stimulates aldosterone release.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Angiotensinogen expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .