

产品名称: Annexin A2 (5A11) Rabbit Monoclonal Antibody
产品货号: AMRe06916



产品概述 (Summary)

产品名称 (Production Name)	Annexin A2 (5A11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA2
别名 (Alternative Names)	Annexin A2; Annexin II; Annexin-2; Calpactin I heavy chain; PAP-IV; Protein I; p36; ANXA2; ANX2; ANX2L4; CAL1H; LPC2D; Annexin2; HEL S 270;
基因 ID (Gene ID)	302.0
蛋白 ID (SwissProt ID)	P07355.A synthetic peptide of human Annexin A2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:20-1:100,ICC/IF 1:20-1:50,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	39kDa

研究背景 (Background)


产品名称: Annexin A2 (5A11) Rabbit Monoclonal Antibody
产品货号: AMRe06916

Originally identified as a protein inhibitor of phospholipase A2, annexin A2 has subsequently been shown to interact with an array of protein and non-protein partners, including F-actin, spectrin, SNARE complexes, RNA, and virus particles. Calcium-regulated membrane-binding protein whose affinity for calcium is greatly enhanced by anionic phospholipids. It binds two calcium ions with high affinity. May be involved in heat-stress response. Inhibits PCSK9-enhanced LDLR degradation, probably reduces PCSK9 protein levels via a translational mechanism but also competes with LDLR for binding with PCSK9 (PubMed:18799458, PubMed:24808179, PubMed:22848640).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of Annexin A2 in C6,3T3,HeLa cell lysates using Annexin A2 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .