

产品名称: Annexin A7 (11B18) Rabbit Monoclonal Antibody
产品货号: AMRe06917



产品概述 (Summary)

产品名称 (Production Name)	Annexin A7 (11B18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA7
别名 (Alternative Names)	Annexin A7; Annexin VII; ANX7; ANXA7; SNX; Synexin;
基因 ID (Gene ID)	310.0
蛋白 ID (SwissProt ID)	P20073.Recombinant protein of human Annexin VII

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,IP 1:50-1:100,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	53kDa

研究背景 (Background)

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

产品名称: Annexin A7 (11B18) Rabbit Monoclonal Antibody
产品货号: AMRe06917




Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of Annexin VII in K562,Rat Brain,HeLa cell lysates using Annexin VII antibody(1:1000 diluted).

注意事项 (Note)

For research use only .