

产品名称: APLP2 (17Z11) Rabbit Monoclonal Antibody
产品货号: AMRe07014



产品概述 (Summary)

产品名称 (Production Name)	APLP2 (17Z11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APLP2
别名 (Alternative Names)	Aplp2; APPH; APPL2; CDEBP;
基因 ID (Gene ID)	334.0
蛋白 ID (SwissProt ID)	Q06481.A synthetic peptide of human APLP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:200-1:500,IF-P 1:200-1:500
蛋白分子量 (Molecular Weight)	87kDa

研究背景 (Background)

May play a role in the regulation of hemostasis. The soluble form may have inhibitory properties towards coagulation


产品名称: APLP2 (17Z11) Rabbit Monoclonal Antibody
产品货号: AMRe07014



factors. May interact with cellular G-protein signaling pathways. May bind to the DNA 5'-GTCACATG-3'(CDEI box). May play a role in the regulation of hemostasis. The soluble form may have inhibitory properties towards coagulation factors. May interact with cellular G-protein signaling pathways. May bind to the DNA 5'-GTCACATG-3'(CDEI box). Inhibits trypsin, chymotrypsin, plasmin, factor XIA and plasma and glandular kallikrein. Modulates the Cu/Zn nitric oxide-catalyzed autodegradation of GPC1 heparan sulfate side chains in fibroblasts (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of APLP2 in HeLa cell lysates using APLP2 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .