

产品名称: ApoA (4M11) Rabbit Monoclonal Antibody
产品货号: AMRe07017



产品概述 (Summary)

产品名称 (Production Name)	ApoA (4M11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LPA
别名 (Alternative Names)	AK38; APOA; apolipoprotein(a); LP; lp(a); LPA;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P08519.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000
蛋白分子量 (Molecular Weight)	501kDa

研究背景 (Background)


产品名称: ApoA (4M11) Rabbit Monoclonal Antibody
产品货号: AMRe07017



Lipoprotein(a) (Lp(a)) is a lipoprotein subclass assembled in the blood from low density lipoprotein (LDL) molecules and apolipoprotein-a (apo-a). Lp(a) recruits inflammatory cells through interaction with Mac-1 integrin. Apo(a) is the main constituent of lipoprotein(a) (Lp(a)). It has serine proteinase activity and is able of autoproteolysis. Inhibits tissue-type plasminogen activator 1. Lp(a) may be a ligand for megalin/Gp 330.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ApoA expression in Human plasma lysate.

注意事项 (Note)

For research use only .